



## Application Areas

- Laboratories: laboratory material, precision instruments, pipettes, sieves, etc.
- Medicine: instrumentation in general, surgical material, etc.
- Odontology: dental prosthesis, instrumentation, etc.
- Optics: instrumentation, optical, frames, etc.
- Industry: printed circuits, Electronic components, etc.
- Jewellery: watches, jewels, etc.
- Automation: Injector clearing.
- Degasification or liquid dissolution.
- Compact substances degermation

## Specifications

Model	UE10SFD
Capacity(L)	10
Internal Size (L x W x H) mm	300×240×150
Overall Size (L x W x H) mm	360×270×310
Working Power (W)	250
Frequency (Khz)	40
Timer (min)	1~99min
Drain	Yes

### WIGGENS GmbH

Gässlesweg 22-24, 75334 Straubenhardt, Germany  
Tel.: 0049 7248 4529088

### WIGGENS China

Room 303, Hall C, Office Building M8, No.1 Jiuxianqiao East  
Road, Chaoyang District, Beijing 100015, China  
Tel: +86 400-809-2068  
Fax: +86 400-809-2068-112  
Email: info@wiggens.com service@wiggens.com  
Website: www.wiggens.com